

Lumify v1.8 Device Compatibility List



Supported Devices*

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy Tab S4 (Model SM-T830)	10.5 in	Wi-Fi	64 GB	5 h 30 min	C	Yes
Samsung Galaxy Tab A -2018 (Model SM-T590)	10.5 in	Wi-Fi	32 GB	3 h 38 min	C	Yes
Samsung Galaxy S9 (SM-G960U1)	5.8 in	Wi-Fi / LTE	32 GB	3 h 45 min	C	Yes
Samsung Galaxy S9+ (SM-G965U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Samsung Galaxy Note9 (SM-N960U1)	6.25 in	Wi-Fi / LTE	128 GB	3 h 26 min	C	Yes
Samsung Galaxy Note8 (SM-N950U1)	6.3 in	Wi-Fi / LTE	64 GB	3 h 30 min	C	Yes
Samsung Galaxy S8 (SM-G950U1)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy S8+ (SM-G955U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Google Pixel 3 / Pixel 3 XL	5.5 in/ 6.3 in	Wi-Fi	64 GB	3 h 43 min	C	Yes
Google Pixel 2 XL	6.0 in	Wi-Fi/LTE	64 GB	3 h	C	Yes
Huawei Mediapad M5 (Model SHT-W09, SHT-AL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	64 GB	3 h 30 min	C	Yes
Xiaomi MiPad 4 (Model M1806D9W) – CHINA ONLY	8.0 in	Wi-Fi / LTE	64 GB	4 h 35 min	C	Yes
Lenovo Smart Tab P10 (Model TB-X705F)	10.1 in	Wi-Fi	32 GB	4 h 53 min	C	Yes
Lenovo Smart Tab M10 (Model TB-605F)	10.1 in	Wi-Fi	32 GB	3 h 43 min	C	Yes
Nokia 6 (Model TA-1003)	5.5 in	Wi-Fi/LTE	32 GB	2 h 30 min	B	Yes
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	32 GB	3 h 30 min	B	Yes
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	Wi-Fi	32 GB	3 h	B	Yes
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.
239034 Rev S March 2019